

Preparation of Chemically Competent *E. coli* -- Rubidium Chloride Method

- Transfer to centrifuge tubes and ice for 20-30 min ↓
1. Inoculate a single colony from a rich plate (Luria-Bertani agar) into 5 ml of rich broth (Luria-Bertani; LB) in a test tube. Shake overnight at 37°C.
 2. Subculture the overnight 1:100 in LB+20 mM MgSO₄ (typically 2.5 ml culture into 250 ml of LB + 5 ml 1M MgSO₄). Grow to OD₅₉₀=0.4-0.6 or Klett=60 (~2-3 h). ~~40-50-0.25-0.4 in plate reader - 200ul sample~~
 3. Centrifuge ^{3220 g} 4,000-rpm 10 min at 4°C. Use 50ml conical tubes. Divide culture evenly into 8 tubes. Eppendorf Centrifuge 5810R - 15 amp version. Swinging bucket rotor.
 4. Gently resuspend pellet in 1/2.5 Volume Unit ice cold TFBI. For 250 ml subculture, use 100 ml TFBI.
 1. discard supernatant ³⁰⁰
 2. add 25 ml TFBI to 4 of the tubes ¹²⁰
 3. resuspend using the 25ml autopipetter
 4. transfer cells to tubes without TFBI and resuspend
 You now have 4 tubes of cells.

Keep all steps on ice and chill all pipets, tubes, flasks, etc. from this point on.

5. Incubate on ice for 5 min.
6. Centrifuge ^{3220 g} 4,000-rpm 10 min 4°C
7. Resuspend pellet in 1/25 original volume cold TFB2. For 250 ml of original subculture, use 10 ml TFB2.
 1. discard supernatant ³⁰⁰
 2. add 10 ml TFB2 to one of the tubes ¹²
 3. resuspend
 4. transfer cells to next tube and resuspend - repeat this step until all 4 tubes have been resuspended
8. Incubate on ice 15-60 min. before aliquoting for storage at -70°C.
9. While incubating, get a tray of ice, label the microfuge tubes (top), and let them chill.
10. Aliquot 110ul and 220ul - about half and half.
11. Quick-freeze the tubes using an EtOH/dry ice bath.
12. Store in -80C freezer.

After the tubes are well-frozen they can be dumped loose into a box or ice-cream carton, or transferred to slots in a storage box. Be careful not to get alcohol on the lips of the tubes. Liquid nitrogen can also be used.

This procedure works with most strains and should routinely give > 10⁷ cfu/ug of pBR322 with reasonably healthy K-12 derivatives (using 0.1 ng/transformation). Frozen cells last at least a year.

Recipes

RB (Luria-Bertani medium)

per liter:

10 g Tryptone (Difco)
5 g Yeast Extract (Difco)
5 g NaCl
2 ml 1N NaOH

TFBI

1.2g 30 mM KOAc (potassium acetate)
4.95g 100 mM RbCl
0.455g 10 mM CaCl₂
2.58g 50 mM MnCl₂
15% glycerol

$$\begin{aligned} 30 \text{ mmol} \times \frac{1 \text{ mol}}{1000 \text{ mmol}} \times 0.4 \text{ L} \times 98.14 \text{ g/mol} &= 1.2 \text{ g KOAc} \\ 100 \text{ mmol} \times \frac{1 \text{ mol}}{1000 \text{ mmol}} \times 6.4 \text{ L} \times 120.92 \text{ g/mol} &= 4.95 \text{ g RbCl} \\ 10 \text{ mM CaCl}_2 \times \frac{1}{1000} \times 0.4 \text{ L} \times 110.99 \text{ g/mol} &= 0.455 \text{ g CaCl}_2 \\ 50 \text{ mM MnCl}_2 \times \frac{1}{1000} \times 0.4 \text{ L} \times 125.84 \text{ g/mol} &= 2.58 \text{ g MnCl}_2 \end{aligned}$$

250ml Culture

0.294 g KOAc Potassium Acetate
1.2092 g RbCl
0.111 g CaCl₂
0.999 g MnCl₂ · 4H₂O
15 ml Glycerol

100ml culture

0.72 g KOAc
2.97 g RbCl
0.273 g CaCl₂
1.548 g MnCl₂

↓
60 ml → 9 ml Glycerol

Adjust to pH 5.8 with acetic acid and filter (0.45 um, Nalgene units or Millipore filters) to sterilize.

It is convenient to make this as:

5 g	RbCl (Alfa)
12.3 ml	KOAc 1 M
4.1 ml	CaCl ₂ 1 M
20.5 ml	MnCl ₂ 1 M (this is pink)
61.5 ml	glycerol
pH to 5.8 with	8 ml HOAc 0.1 M

~~1.25g RbCl~~

make up to 410 ml; distribute in 100 ml sterile aliquots; and use 1 aliquot/250 ml culture.

TFBII

1.25g 10 mM MOPS or PIPES
0.18g 75 mM CaCl₂
10 mM RbCl
15% glycerol

$$\begin{aligned} 10 \times \frac{1}{1000} \times 0.15 \times 209 &= 0.31 \text{ g MOPS} \\ 75 \text{ mmol} \times \frac{1 \text{ mol}}{1000 \text{ mmol}} \times 10.150 \text{ L} \times 110.99 \text{ g/mol} &= 1.25 \text{ g CaCl}_2 \\ 10 \text{ mM RbCl} \times \frac{1}{1000} \times 0.15 \times 120.92 \text{ g/mol} &= 0.18 \text{ g RbCl} \end{aligned}$$

Adjust pH to 6.5 with KOH and filter to sterilize

Make up as

1.5ml	MOPS 1 M pH 6.5 (this is yellow)
11.25ml	CaCl ₂ 1 M
1.5ml	RbCl 1M
22.5 ml	glycerol
pH with	1 N KOH

make to 150 ml; filter; use 10 ml per original 250 ml culture.

Procedure adapted from the John Innes Institute (Norwich, England) via Joseph Utermohlen (Univ. of Arizona).

This procedure improves transformation with most strain backgrounds, when compared with CaCl₂ procedure.